
PNEUMASOX[®]

COMPRESSION THERAPY

INSTRUCTIONS FOR USE

A NEW WAY OF PERFECT COMPRESSION

PNEUMASOX

COMPRESSION THERAPY

INDEX

NL	Gebruiksaanwijzing	4 - 7
EN	Instructions for use	8 - 11
DE	Gebrauchsanweisung	12 - 15
FR	Mode d'emploi	16 - 19
IT	Istruzioni per l'uso	20 - 23
ES	Manual de instrucciones	24 - 27
PT	Manual do utilizador	28 - 31
SK	Návod na použitie	32 - 35

BENAMING

PneumaSox

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PneumaSox

PRODUCT DESCRIPTION

PneumaSox is a static pneumatic compression device for compression therapy, to be used on the lower leg instead of bandages. It utilises air pressure. The pressure exerted on the leg by the cuff prevents oedema formation, promotes venous return and supports lymphatic drainage.

COMPOSITION

PneumaSox has two components:

A calf-length compression cuff

- This consists of a flat, air chamber that can be wrapped around the leg and secured with Velcro.
- The cuff is made of Polyamide / Polyurethane elastomer.

And a hand-operated pump with pressure sensor.**CHARACTERISTICS**

- PneumaSox is a static pneumatic compression device for compression therapy, to be used on the lower leg instead of bandages.
- PneumaSox uses a pressure sensor to provide precise, optimal pressure, evenly distributed over the leg.
- This pressure can easily be set, and adjusted.
- The even distribution of pressure contributes to a high degree of comfort when wearing the cuff.
- PneumaSox can be worn with shoes (see figure 6).

INDICATIONS/APPLICATION AREAS

PneumaSox compression therapy is designed to provide static, pneumatic compression to the lower leg, for treatment of:

- Ulcus cruris venosum (end-stage CVI)
- Chronic venous insufficiency (CVI)
- Superficial and deep venous thrombosis
- Dependency oedema
- Erysipelas
- Arterial ulcers (after pressure measurement)
- Martorell ulcer
- Traumatic ulcers
- Diabetic ulcers
- Other causes of ulcers associated with oedema.

CONTRAINDICATIONS

Absolute contraindications for the application of compression therapy:

- Arterial insufficiency (distal systolic pressure (DSP) < 60mmHg and/or ankle/arm index < 0.5) and total occlusion of the deep venous system
- Severe heart failure (NYHA class IV)
- Serious skin abnormalities
- Allergy regarding one of the components of the compression material
- When the patient is incapable of operating PneumaSox.

Figure 1

Figure 2

Figure 3

BEFORE USE

- Always follow the advice of your healthcare professional for using PneumaSox. Their advice takes precedence.
- In general, it is recommended that you use PneumaSox from the time you get up until you go to bed.

DONNING THE CUFF

- Make sure the cuff is empty when you put it on.
- The cuff **labelled L** is for the left leg and the cuff **labelled R** is for the right leg.
- The cuff must not come into contact with damaged skin.
You can wear an understocking under the cuff, in consultation with your doctor. In such case, make sure that such understocking has a smooth, wrinkle-free fit. If an understocking worn under the cuff is too long, fold it down over the cuff, after putting on the cuff.
- Insert your foot into the cuff (**see figure 1**).
- Close the cuff around the leg. Make sure the Velcro straps are properly closed (**see figure 2**).
- After putting on the cuff, the (Velcro) closure should be on the outside of the leg.

INFLATING THE CUFF / USING THE PRESSURE SENSOR

- Connect the hand pump with pressure sensor to the valve of the cuff. Once the cuff has been put on, this valve can be found on the inside of your leg (**see figure 3**).
- Switch on the pressure sensor by pressing the button for 1 second. Two green lights will start flashing.
- Increase the pressure in the cuff using the hand pump (**see figure 4**). 'Pressure' means how hard the cuff presses against your leg.
- The green lights will continue to flash while the pressure is too low. Keep pumping air into the cuff, for as long as the green lights flash. The green lights will flash faster as the pressure increases. Once the cuff reaches the correct pressure, the lights will stop flashing and remain continuously on (**see figure 5**). Stop pumping at that time.
- A little excess pressure, that might result from one extra squeeze of the hand pump, does no harm. In such case, the red lights will then flash. If the pressure is unacceptably high, the red lights will stop flashing and stay on continuously. In such case, you must reduce the pressure in the cuff. This can be done by moving your toes up and down (*see: After inflation*). When the pressure is right, the green lights will stay on continuously.
- If the red lights remain on continuously, let some air out of the cuff (*see: Taking of the cuff*), then pump air in again until the green lights are continuously on. The pressure sensor switches itself off after use.

AFTER INFLATION

- Once the cuff is inflated, move your toes up and down 10 times. You will notice that this slightly reduces the pressure in the cuff.
- After this, re-inflate the cuff (*see: Inflating the cuff/ Using the pressure sensor*).
- After inflating, carefully disconnect the pump with pressure sensor from the PneumaSox. Do this by grabbing the end of the tube (where it is attached to the PneumaSox) and pulling gently.

Figure 4

Figure 5

Figure 6

CHECKING THE PRESSURE

- An increase or decrease in the swelling of the leg may increase or decrease the pressure exerted by the cuff. It is therefore important to check the pressure regularly, throughout the day, and adjust it if necessary.
- To do this, connect the hand pump with pressure sensor to the valve on the cuff (see **Figure 3**) and switch on the pressure sensor.
- Check the pressure 1 hour after putting on the cuff, around lunchtime and before dinner. Checking the pressure more often than this is recommended.

TAKING OFF THE CUFF

- Open the Velcro straps and remove the cuff from your leg.
- Expel as much air as possible from the cuff. To do this, insert the transparent nozzle at the back of the hand pump into the valve of the cuff. Now gently roll up the cuff, using a little pressure, to expel the air from the cuff.

WEARING SHOES

- Put on your shoes after putting on the cuff.
- Make sure that the lower end of the cuff, at your ankle, is properly seated in your shoe.
- Then inflate the cuff to the set pressure (see: *Inflating the cuff/ using the pressure sensor* + **figure 4**).
- After removing your shoes, check the pressure in the cuff and pump it up again if necessary (see: *Checking the pressure*).
- If the cuff is partially or fully inflated, it can be difficult to put on your shoes. If so, deflate the cuff slightly. After putting on your shoes, you can re-inflate the cuff.

WASHING INSTRUCTIONS

Cuff

- The cuff can be machine washed at 40° centigrade.
- The cuff may also be washed by hand (rinse well and do not wring out).
- The Velcro straps must be closed before washing!
- Wash cuff separately from other laundry items the first time you wash it, as it may release colour.
- Placing the cuff in a washing net provides extra protection.
- Do not use fabric softener or bleach.
- The cuff must not be chemically or dry-cleaned.

Hand pump

- If necessary, the hand pump can be cleaned with a damp cloth, using a mild dishwashing detergent.
- Do not immerse the pump in water or liquid - this will damage the pump.

DRYING THE CUFF

- The cuff will dry quickly. Hang the cuff up to dry. Do not stretch it or peg it.
- Do not dry the cuff in a tumble dryer, on a radiator, or in the sun.
- Do not iron or steam the cuff.

SETTING THE PRESSURE SENSOR

- **THE CORRECT PRESSURE SHOULD ONLY BE SET BY A HEALTHCARE PROFESSIONAL.**
- The pressure sensor is pre-set by the factory to a preferred pressure of 40 mmHg.
- Once this 'preferred pressure' is reached, the green lights will stop flashing and remain on continuously. The preferred pressure can easily be adjusted before use, by a healthcare professional. Adjustment is possible in increments of 10 mmHg to 100 mmHg.
- Switch on the sensor. To set the pressure, press the button continuously for 7 seconds. The green lights will illuminate brightly for a moment.
- Now press the button as many times as the desired pressure divided by 10. For example, 6 times for a desired pressure of 60 mmHg. Wait a few seconds.
- The lights will flash a number of times to enable checking of the desired pressure (i.e. they will flash n times where $n \times 10 =$ the desired, set pressure in mmHg). The sensor will thereafter enter the operating mode.

FOR HOW LONG SHOULD I CONTINUE TO USE PNEUMASOX?

- Compression therapy using PneumaSox is usually long-term.
- Consult your doctor on the duration of the therapy.

USAGE PERIOD

The useful life of PneumaSox is 6 months. Useful life depends on proper handling of the product (e.g. how it is maintained, put on and taken off) and the number of times the product is used.

SPECIAL ATTENTION AND PRECAUTIONS

- Keep out of reach of children.
- Use only on the advice of your doctor, nurse or home carer.
- The pre-set, preferred pressure for PneumaSox may only be set by your doctor, nurse or home carer.
- The cuff should never hurt you or be uncomfortable to wear. You should feel a firm but comfortable pressure. If you experience severe discomfort, let some air out of the compression part. If this does not improve your comfort, remove the cuff and consult your doctor.
- Do not use PneumaSox if you have poor, obstructed circulation. In such case, consult your doctor first. Compression applied to a leg that already suffers from partially obstructed blood supply can lead to potentially serious complications.
- Each PneumaSox cuff can only be worn by one user, for hygienic reasons.
- Exercise care if ambient air pressure changes (such as in an airplane or in the mountains), as this will change the compression exerted by the cuff. If you feel a noticeable change in compression, release air from the cuff and re-inflate it to the set pressure.
- Follow all wound care instructions issued by your healthcare professional.
- Please note: in case of damaged skin, always use an understocking under the cuff.
- Do not use in an MRI environment or hyperbaric chamber.
- Only use the pressure sensor and hand pump supplied in the package containing the PneumaSox.
- Do not swim or bathe while wearing PneumaSox.
- Do not use if the packaging is damaged or already open on receipt of the product, unless in consultation with your healthcare professional.

SIDE EFFECTS

- When used correctly, there are no known side effects.
- However, if a negative change (e.g. skin irritation or worsening of your symptoms) occurs while using PneumaSox, stop using the product and contact your doctor/healthcare professional immediately.
- Wearing the cuff in direct contact with your skin may cause irritation. In such case, wear an understocking under the cuff.
- If your leg hurts or your foot becomes pale, cold or numb, remove the cuff from your leg.
- Serious events related to this device should be reported to the manufacturer and to the competent authorities (of the Member State) where the user is domiciled.

STORAGE

- Product and storage place must be dry.
- PneumaSox can be stored for up to 60 months. This storage period relates to the battery life of the pressure sensor.

DISPOSAL

- The cuff may be disposed of via the usual waste disposal process. There are no special criteria for its disposal as waste.
- The hand pump with pressure sensor contains a battery and must be disposed of as household hazardous waste.

DATE OF LATEST REVISION OF THESE INSTRUCTIONS FOR USE

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BAP Medical B.V.