



Setting the course for better healing

**Personalized Transient Skin-like Layer
Anywhere · Anytime**

The Spincare is a portable, handheld, easy to use electrospinning device for the creation of a fully personalized, in-situ transient skin-like layer for the treatment of wounds and burns.

The Spincare System consists of a handheld device and a sterile pre-filled disposable ampule which creates a customized nano-fibrous transient skin-like layer that remains on the wound until the skin underneath is fully epithelialized.

Spincare is the first and only commercialized available portable system that integrates electrospinning technology offering immediate wound care treatment that outperforms traditional methods.




Highlights of the Spincare System

Live While You Heal



No contact
one-time
application




Allows
early showers



Nano-fibrous,
multi-layer, transparent
transient skin-like layer




Pain free excellent
adherence to
wound bed




Excellent healing
and scarring



Mobility is
not impaired



No painful
redressing



Cost effective,
reduced financial burden



Personalized to fit
all shapes including
hard to dress areas

About Nanomedic

Nanomedic Technologies Ltd. is an Israeli bio-medical high-technology company specializing in the development and marketing of its Spincare portable technology; a patented platform technology with wide market applications both within and outside the hospital environments.

Nanomedic is focused on the development of next generation wound care solutions and technologies for the benefit of patients, medical staff and the healthcare system.

Contact

Nanomedic Technologies Ltd.

Global Park, 2 Yodfat St.
Northern Ind. Zone, Lod
Israel 7129104

t +972 (0)8-915-3001

e info@nanomedic.com

nanomedic.com

